

# **FastGene**<sup>®</sup> Scriptase II

## cDNA Synthesis 5x ReadyMix

Manual

**LS64**

Version 1.0

April 2018

*For Research use only*

Store at -20°C



**NIPPON Genetics EUROPE**

 [www.nippongenetics.eu](http://www.nippongenetics.eu)  
 [info@nippongenetics.eu](mailto:info@nippongenetics.eu)

**TEL** +49 2421 55 4960  
**FAX** +49 2421 55 49611

## Product description

### Reverse transcription

The ability to reverse the process of translation from DNA to RNA has been a crucial application of modern molecular biology. The enzyme, a RNA-dependent DNA-Polymerase is called reverse transcriptase. Reversing the messenger RNA delivers the complementary DNA, which can be detected or quantified using PCR and qPCR approaches.

The lack of the 3'→5' exonuclease activity made the MuLV increasingly the enzyme of choice. The MuLV has also been used as a template to engineer new enzymes with additional features.

The FastGene® Scriptase II 5x ReadyMix is an engineered enzyme with a reduced RNase H activity and comes ready-to-use with all necessary ingredients to perform a reverse transcription.

## Composition (100 Reactions)

- FastGene® Scriptase II
- 5x FastGene® Scriptase II buffer
- dNTP Mixture
- RNase Inhibitor
- Random Hexamer
- Helper Protein

## Quality Control

- Purity: >99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free
- Inhibitor-free
- Satisfactory Yield of DNA products

# Protocol

1. Mix the ReadyMix with the RNA Template using the following setup:

5x FastGene® Scriptase II ReadyMix (LS64)		4 $\mu$ l
High quality RNA-template*		x $\mu$ l
Sterile Water	up to	20 $\mu$ l

\*Concentration of template RNA can be up to 2.5  $\mu$ g

2. Incubate at 25 °C for 10 minutes.

3. Incubate at 42 °C for 60 minutes\*.

**Note:** The recommended incubation time can be reduced to 5 minutes depending on the product abundance and size. Please download our Technical Note showing the performance at lower temperature (down to 5 min.) at our website [www.nippongenetics.eu/en/](http://www.nippongenetics.eu/en/).



3. Inactivate the FastGene Scriptase II by incubating at 85°C for 5 min.

4. Store the cDNA at -20 °C or use it for downstream application.

**Note:** Up to 10 % of the qPCR reaction volume can be of the undiluted synthesized cDNA .

# **FastGene**® Scriptase II

## cDNA Synthesis 5x ReadyMix

Please contact us for additional information:

**info@nippongenetics.eu**

Please contact our technical support:

**support@nippongenetics.eu**

LS52	<b>FastGene® Scriptase Basic</b> (20.000 units at 200 U/μl)	100 Reactions
LS62	<b>FastGene® Scriptase Basic cDNA Synthesis</b> containing Oligo dTs, random hexamer and RNase inhibitor	100 Reactions
LS53	<b>FastGene® Scriptase II</b> (20.000 units at 200 U/μl)	100 Reactions
LS63	<b>FastGene® Scriptase II cDNA Synthesis</b> containing Oligo dTs, random hexamer and RNase inhibitor	100 Reactions
LS64	<b>FastGene® Scriptase II cDNA Synthesis 5x ReadyMix</b> containing random hexamer, RNase inhibitor and Helper Protein	100 Reactions

 **NIPPON Genetics EUROPE**

 [www.nippongenetics.eu](http://www.nippongenetics.eu)  
 [info@nippongenetics.eu](mailto:info@nippongenetics.eu)

**TEL** +49 2421 55 4960  
**FAX** +49 2421 55 49611